



D-[1-¹³C]allose
¹³C {¹H} Spectrum (75 MHz) in D₂O



Assignments

α -pyranose
C1 94.3
C2 68.6
C3 73.2
C4 67.6
C5 68.3
C6 62.3

β -pyranose
C1 94.9
C2 72.8
C3 72.7
C4 68.3
C5 75.1
C6 62.8

α -furanose
C1 97.5

β -furanose
C1 102.3

Specialty Stable Isotopes for the Biomedical Sciences

115 South Hill Street
South Bend, Indiana, USA 46617-2701
www.omicronbio.com

Phone: 574-287-6910
Fax: 574-287-7165
Email: omicron@omicronbio.com