


**D-galactose**  
 **$^{13}\text{C}$  { $^1\text{H}$ } Spectrum (75 MHz) in  $\text{D}_2\text{O}$**



**Assignments**

$\alpha$ -pyranose  
C1 93.6  
C2 69.8  
C3 70.7  
C4 70.8  
C5 72.0  
C6 62.7

$\beta$ -pyranose  
C1 97.9  
C2 73.4  
C3 74.3  
C4 70.2  
C5 76.6  
C6 62.5

$\alpha$ -furanose  
C1 96.5

$\beta$ -furanose  
C1 102.5

---

**Specialty Stable Isotopes for the Biomedical Sciences**

115 South Hill Street  
South Bend, Indiana, USA 46617-2701  
www.omicronbio.com

Phone: 574-287-6910  
Fax: 574-287-7165  
Email: omicron@omicronbio.com