


D-[1-¹³C]gulose
¹³C {¹H} Spectrum (75 MHz) in D₂O




Assignments

α -pyranose

C1	94.4
C2	66.3
C3	72.4
C4	71.0
C5	68.0
C6	62.5

β -pyranose

C1	95.4
C2	70.7
C3	72.8
C4	71.0
C5	75.3
C6	62.6



Specialty Stable Isotopes for the Biomedical Sciences

115 South Hill Street
South Bend, Indiana, USA 46617-2701
www.omicronbio.com

Phone: 574-287-6910
Fax: 574-287-7165
Email: omicron@omicronbio.com