



**D-[1-<sup>13</sup>C]mannose**  
**<sup>13</sup>C {<sup>1</sup>H} Spectrum (75 MHz) in D<sub>2</sub>O**



**Assignments**

$\alpha$ -pyranose	
C1	95.5
C2	72.2
C3	71.7
C4	68.4
C5	73.9
C6	62.5

$\beta$ -pyranose	
C1	95.2
C2	72.7
C3	74.5
C4	68.1
C5	77.6
C6	62.5



**Specialty Stable Isotopes for the Biomedical Sciences**

115 South Hill Street  
South Bend, Indiana, USA 46617-2701  
www.omicronbio.com

Phone: 574-287-6910  
Fax: 574-287-7165  
Email: omicron@omicronbio.com