



D-[1-¹³C]talose
¹³C {¹H} Spectrum (75 MHz) in D₂O



Assignments

α -pyranose

C1	96.1
C2	72.2*
C3	71.1*
C4	66.6
C5	72.6
C6	63.0

β -pyranose

C1	95.6
C2	73.0*
C3	70.1*
C4	69.9
C5	77.1
C6	62.7

α -furanose

C1	102.4
C2	76.7
C3	73.3
C4	83.3
C5	72.1
C6	64.3

β -furanose

C1	98.0
C2	72.1
C3	72.6
C4	83.9
C5	72.3
C6	64.4

Specialty Stable Isotopes for the Biomedical Sciences

115 South Hill Street
South Bend, Indiana, USA 46617-2701
www.omicronbio.com

Phone: 574-287-6910
Fax: 574-287-7165
Email: omicron@omicronbio.com