


D-[1-¹³C]xylose
¹³C {¹H} Spectrum (75 MHz) in D₂O




Assignments

α-pyranose
C1 93.7
C2 73.0
C3 74.3
C4 70.9
C5 62.4

β-pyranose
C1 98.1
C2 75.5
C3 77.3
C4 70.7
C5 66.7

α-furanose
C1 97.0

β-furanose
C1 103.2



Specialty Stable Isotopes for the Biomedical Sciences

115 South Hill Street
South Bend, Indiana, USA 46617-2701
www.omicronbio.com

Phone: 574-287-6910
Fax: 574-287-7165
Email: omicron@omicronbio.com