


D-lyxose
 ^{13}C { ^1H } Spectrum (75 MHz) in D_2O



Assignments

α -pyranose
C1 95.5
C2 71.5
C3 72.0
C4 69.0
C5 64.6

β -pyranose
C1 95.9
C2 71.5
C3 74.2
C4 68.0
C5 65.7

Specialty Stable Isotopes for the Biomedical Sciences

115 South Hill Street
South Bend, Indiana, USA 46617-2701
www.omicronbio.com

Phone: 574-287-6910
Fax: 574-287-7165
Email: omicron@omicronbio.com