


D-[1-¹³C]glucose
¹³C {¹H} Spectrum (75 MHz) in D₂O





Assignments

α -pyranose

C1	93.6
C2	73.0
C3	74.2
C4	71.14
C5	72.9
C6	62.1

β -pyranose

C1	97.4
C2	75.6
C3	77.2
C4	71.09
C5	77.4
C6	62.3



Specialty Stable Isotopes for the Biomedical Sciences

115 South Hill Street
South Bend, Indiana, USA 46617-2701
www.omicronbio.com

Phone: 574-287-6910
Fax: 574-287-7165
Email: omicron@omicronbio.com